

Metabolomic Data Analysis with MetaboAnalystR

Name: guest623839008349489450

February 7, 2018

1 Background

Understanding the functional importance of metabolites in untargeted metabolomics is limited due to challenges with metabolite identification. To reduce problems associated with compound misidentification and thereby pathway misinterpretation is to shift the unit of analysis from individual compounds to individual pathways. In particular, the mummichog algorithm (5) offers an elegant and efficient implementation of this concept. Mummichog bypasses the bottleneck of metabolite identification prior to pathway analysis by leveraging a priori pathway and network knowledge to directly infer biological activity based on MS peaks. Due to its popularity and repeated user requests, we have implemented the mummichog algorithm (version 1.0.10) from Li et al. 2013, which has been carefully translated from the Python programming language to R, and includes an expanded knowledgebase of 21 organisms for pathway analysis. In particular, this module by-passes the bottle-neck of metabolite identification prior to pathway analysis, leveraging a priori knowledge from genome-scale metabolic models and KEGG metabolic pathways. For instance, conventional approaches require statistically significant metabolites to be identified prior to pathway analysis, which can be incredibly time-consuming, whereas here, m/z features are used in conjunction with metabolic models/pathways to directly infer pathway level knowledge from m/z features. For further details, please refer to Li et al. 2013 (PMC3701697).

2 Overview

The MS Peaks to Pathways module consists of three steps - uploading the user's data, selection of a pathway library, and pathway analysis.

3 Data Input

The MS Peaks to Pathways module accepts either a list of significant m/z features, a list of all m/z features, or a peak intensity table. The format of the data must be specified, identifying whether the samples are in rows or columns, and whether or not the data is paired. The data may either be .csv or .txt files.

3.0.1 MS Peaks to Pathways: Reading Data

The data must be uploaded as a three column table containing the m/z features, p-values, and statistical scores (e.g. t-scores or fold-change values).

A total of 3934 m/z features were found in your uploaded data.

3.0.2 Parameters

Users also need to specify the mass accuracy, the ion mode (positive or negative), and the p-value cutoff to delineate between significantly enriched and non-significantly enriched m/z features. Currently, Metabo-Analyst 4.0 only supports the handling of peaks obtained from high-resolution MS instruments such as Orbitrap, or Fourier Transform (FT)-MS instruments as recommended by the original mummichog implementation.

The selected mass accuracy of your MS instrument in ppm is: ten; The selected mode of your MS instrument is: positive; The selected p-value cutoff is: 1e-04.

3.0.3 Library

The knowledge-base for this module consists of five genome-scale metabolic models obtained from the original Python implementation which have either been manually curated or downloaded from BioCyc, as well as an expanded library of 21 organisms derived from KEGG metabolic pathways. Users must select one of 21 KEGG pathway libraries, or one of five metabolic models.

The user's selected library is: *hsa_mfn*.

4 Output

The aim of this module is to leverage the power of known metabolic models/pathways to gain functional insight directly from m/z features. There are three steps in this module, 1) Permutations: A list of metabolites (the same length as the number of significant m/z features) are inferred from the user's uploaded set of m/z features, considering all potential matches (isotopes/adducts). These tentative compounds are then mapped onto known metabolic pathways for the selected organism. For each pathway, a fisher's exact or hypergeometric p-value is calculated. 2) Step 1 is repeated multiple times to calculate the null distribution of p-values for all pathways, and is modeled as a Gamma distribution. 3) Following this, the significant m/z features are used to calculate the p-values for each pathway (Step 1). These p-values are then adjusted for the permutations.

5 Pathway Analysis Results Table

The output of the MS Peaks to Pathways module consists of a table of results containing ranked pathways that are enriched in the user-uploaded data. The table includes the total number of hits, their raw p-values (Fisher's exact test or Hypergeometric), their EASE score, and the p-value modeled on user data using a Gamma distribution.

Table 1: Results of the mummichog pathway analysis

	Pathway total	Hits.total	Hits.sig	EASE	FET	Gamma
Tryptophan metabolism	94.00	64.00	21.00	0.01	0.00	0.00
Ascorbate (Vitamin C) and Aldarate Metabolism	29.00	18.00	9.00	0.01	0.00	0.00
Aminosugars metabolism	69.00	29.00	12.00	0.01	0.00	0.00
Nitrogen metabolism	6.00	4.00	4.00	0.02	0.00	0.01
N-Glycan biosynthesis	48.00	14.00	7.00	0.03	0.01	0.01
Pyrimidine metabolism	70.00	43.00	14.00	0.05	0.02	0.01
Vitamin B3 (nicotinate and nicotinamide) metabolism	28.00	19.00	8.00	0.05	0.02	0.01
Sialic acid metabolism	107.00	28.00	10.00	0.06	0.03	0.01
Alanine and Aspartate Metabolism	30.00	17.00	7.00	0.08	0.03	0.01
Glutathione Metabolism	19.00	10.00	5.00	0.10	0.03	0.01
Hexose phosphorylation	20.00	18.00	7.00	0.10	0.04	0.01
Arginine and Proline Metabolism	45.00	31.00	10.00	0.11	0.05	0.01
Glycosphingolipid biosynthesis - ganglioseries	62.00	11.00	5.00	0.13	0.04	0.01
Glutamate metabolism	15.00	11.00	5.00	0.13	0.04	0.01
Methionine and cysteine metabolism	94.00	46.00	13.00	0.14	0.08	0.01
Purine metabolism	80.00	51.00	14.00	0.14	0.08	0.01
Parathio degradation	6.00	4.00	3.00	0.16	0.02	0.01
Vitamin B9 (folate) metabolism	33.00	12.00	5.00	0.17	0.06	0.01
Glycosphingolipid biosynthesis - globoseries	16.00	8.00	4.00	0.18	0.05	0.01
Starch and Sucrose Metabolism	33.00	15.00	5.00	0.31	0.13	0.01
N-Glycan Degradation	16.00	6.00	3.00	0.31	0.08	0.02
Linoleate metabolism	46.00	21.00	6.00	0.36	0.19	0.02
Urea cycle/amino group metabolism	85.00	42.00	10.00	0.39	0.26	0.02
Glycosphingolipid metabolism	67.00	28.00	7.00	0.43	0.27	0.02
Aspartate and asparagine metabolism	114.00	61.00	13.00	0.48	0.36	0.03
Caffeine metabolism	11.00	9.00	3.00	0.53	0.23	0.04
Lipoate metabolism	8.00	4.00	2.00	0.57	0.16	0.04
Vitamin B1 (thiamin) metabolism	20.00	10.00	3.00	0.59	0.29	0.05
Drug metabolism - cytochrome P450	53.00	38.00	8.00	0.59	0.43	0.05
TCA cycle	31.00	16.00	4.00	0.60	0.36	0.05
Pentose phosphate pathway	37.00	33.00	7.00	0.61	0.43	0.05
Porphyrim metabolism	43.00	22.00	5.00	0.62	0.40	0.05
Alkaloid biosynthesis II	10.00	5.00	2.00	0.65	0.24	0.06
Glycosphingolipid biosynthesis - lactoseries	14.00	5.00	2.00	0.65	0.24	0.06
O-Glycan biosynthesis	16.00	5.00	2.00	0.65	0.24	0.06
Tyrosine metabolism	160.00	87.00	16.00	0.69	0.59	0.07
Glycolysis and Gluconeogenesis	49.00	36.00	7.00	0.70	0.53	0.08
C5-Branched dibasic acid metabolism	10.00	6.00	2.00	0.71	0.32	0.08
Chondroitin sulfate degradation	37.00	6.00	2.00	0.71	0.32	0.08
Blood Group Biosynthesis	44.00	6.00	2.00	0.71	0.32	0.08
Glycosphingolipid biosynthesis - neolactoseries	16.00	6.00	2.00	0.71	0.32	0.08
Keratan sulfate biosynthesis	66.00	6.00	2.00	0.71	0.32	0.08
Glycerophospholipid metabolism	156.00	37.00	7.00	0.72	0.56	0.09
Pentose and Glucuronate Interconversions	15.00	13.00	3.00	0.73	0.46	0.09
Heparan sulfate degradation	34.00	7.00	2.00	0.77	0.39	0.11
Lysine metabolism	52.00	28.00	5.00	0.80	0.63	0.13
Keratan sulfate degradation	68.00	8.00	2.00	0.81	0.46	0.13
Vitamin B5 - CoA biosynthesis from pantothenate	12.00	8.00	2.00	0.81	0.46	0.13
Carbon fixation	10.00	8.00	2.00	0.81	0.46	0.13
Vitamin E metabolism	54.00	16.00	3.00	0.83	0.61	0.15
Drug metabolism - other enzymes	31.00	16.00	3.00	0.83	0.61	0.15
Histidine metabolism	33.00	19.00	3.00	0.90	0.73	0.22
Butanoate metabolism	34.00	27.00	4.00	0.91	0.78	0.24
Xenobiotics metabolism	110.00	70.00	10.00	0.93	0.88	0.29
Squalene and cholesterol biosynthesis	55.00	22.00	3.00	0.94	0.81	0.30
Selenoamino acid metabolism	35.00	22.00	3.00	0.94	0.81	0.30
Glycine, serine, alanine and threonine metabolism	88.00	51.00	7.00	0.94	0.88	0.30
Galactose metabolism	41.00	30.00	4.00	0.94	0.85	0.31
Biopterin metabolism	22.00	14.00	2.00	0.95	0.77	0.32
Saturated fatty acids beta-oxidation	36.00	14.00	2.00	0.95	0.77	0.32
Beta-Alanine metabolism	20.00	14.00	2.00	0.95	0.77	0.32
Putative anti-Inflammatory metabolites formation from EPA	27.00	15.00	2.00	0.96	0.81	0.35
Phosphatidylinositol phosphate metabolism	59.00	24.00	3.00	0.96	0.86	0.35
Carnitine shuttle	72.00	16.00	2.00	0.96	0.84	0.38
Propanoate metabolism	31.00	20.00	2.00	0.98	0.92	0.50
Leukotriene metabolism	92.00	30.00	3.00	0.99	0.94	0.50
Prostaglandin formation from arachidonate	78.00	42.00	4.00	0.99	0.97	0.56
Androgen and estrogen biosynthesis and metabolism	95.00	53.00	5.00	0.99	0.98	0.61
Fructose and mannose metabolism	33.00	26.00	2.00	1.00	0.97	0.63
Valine, leucine and isoleucine degradation	65.00	40.00	3.00	1.00	0.99	0.69
Arachidonic acid metabolism	95.00	44.00	3.00	1.00	0.99	0.74
C21-steroid hormone biosynthesis and metabolism	112.00	57.00	2.00	1.00	1.00	0.93
Hyaluronan Metabolism	8.00	4.00	1.00	1.00	0.57	1.00
Pyruvate Metabolism	20.00	12.00	1.00	1.00	0.92	1.00
Vitamin H (biotin) metabolism	5.00	4.00	1.00	1.00	0.57	1.00
Glyoxylate and Dicarboxylate Metabolism	12.00	6.00	1.00	1.00	0.72	1.00
De novo fatty acid biosynthesis	106.00	19.00	1.00	1.00	0.98	1.00
Fatty Acid Metabolism	63.00	10.00	1.00	1.00	0.88	1.00
CoA Catabolism	7.00	5.00	1.00	1.00	0.65	1.00
D4&E4-neuroprostanes formation	37.00	19.00	1.00	1.00	0.98	1.00
Dynorphin metabolism	8.00	4.00	1.00	1.00	0.57	1.00
Vitamin B6 (pyridoxine) metabolism	11.00	7.00	1.00	1.00	0.77	1.00

6 MS Peaks to Pathway Output: Compound Matching Table

The output of the MS Peaks to Pathways module also consists of a comprehensive table containing the compound matching information for all user-uploaded m/z features. The table has four columns, containing the Query.Mass of each feature, the predicted Matched.Compound for each feature, the Matched.Form, and the Mass.Diff. As the file can be very long (>40 pages), please download it separately on the Downloads page of MetaboAnalyst.

7 Appendix: R Command History

```
[1] "InitDataObjects(\"mass_all\", \"mummichog\", FALSE)"
[2] "mSet<-Read.PeakListData(mSet, \"Replacing_with_your_file_path\");"
[3] "mSet<-UpdateMummichogParameters(mSet, \"ten\", \"positive\", 1.0E-4);"
[4] "mSet<-SanityCheckMummichogData(mSet)"
[5] "mSet<-SetMass.PathLib(mSet, \"hsa_mfn\")"
[6] "mSet<-PerformMummichog(mSet, \"fisher\", \"gamma\")"
[7] "mSet<-SaveTransformedData(mSet)"
```

The report was generated on Wed Feb 7 16:12:35 2018 with R version 3.4.1 (2017-06-30).